

# The Right Sample for Molecular Assays

## Twin Sample Biopsy Needle

*Two sample chambers capture twin samples from one location in one pass.  
No more blind biopsies.*

*Pathology shows what's in the first sample.  
Extract the second sample in one lab and send  
DNA, RNA & protein for analysis —  
No contradictions. Drivers confirmed.*


*The Twin Sample Needle puts pathology and  
molecular analyses on the same page.*

*This directly addresses Recommendation 9  
from the newly released: **Biomarker Tests  
for Molecularly Targeted Therapies:  
Key to Unlocking Precision Medicine**  
from the National Academy of Sciences.<sup>†</sup>*

*“Enhance specimen handling and documentation  
to ensure patient safety and the accuracy of biomarker test results.”*

*“Health care professionals who collect, process, and handle patient biomaterials  
for biomarker testing should ensure that adequate tissue is acquired to perform  
all necessary testing so patients are protected from unnecessary or repeated  
procedures, and that specimens are properly handled...”*

<sup>†</sup> National Academies of Sciences, Engineering, and Medicine. 2016. **Biomarker tests for molecularly targeted therapies: Key to unlocking precision medicine**  
Washington, DC: The National Academies Press. doi: 10.17226/21860



***twinsample.com***  
**AACR 2016 • Booth 2337**

*Epitome Pharmaceuticals Ltd., San Diego, CA. Call toll free (855) 373-8663*